



Nexus on Unicone

Improved cleanability and maintenance




True friction fit



Optional metal framework



Supports a wide variety of restorative options



Available in 3 heights and 2 tapers for optimal fit and positioning

The Nexus on Unicone is a digitally created overdenture structure that brings all the advantages of the Nexus overlays, designed to suit the Unicone conical abutment solution. It is designed using Osteon's proprietary software and precision milled to ensure a true friction fit. These structures are intended to retain a pre-manufactured pure Titanium cap, and supports a wide variety of restorative overlay options.